

Product datasheet for **GA217460**

Mouse Tnip3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnip3 CRISPRa kit - CRISPR gene activation of mouse TNFAIP3 interacting protein 3
Symbol:	Tnip3
Synonyms:	9030611K07Rik; ABIN-3; LIND
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)