

Product datasheet for **GA217481**

Mouse Tmprss9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmprss9 CRISPRa kit - CRISPR gene activation of mouse transmembrane protease, serine 9
Symbol:	Tmprss9
Components:	GA217481G1, Tmprss9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217481G2, Tmprss9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217481G3, Tmprss9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Serase-1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)