


## Product datasheet for **GA217553**

### Mouse Ccdc185 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Ccdc185 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 185                                                                                                                                                                                                                          |
| Symbol:       | Ccdc185                                                                                                                                                                                                                                                                                                          |
| Components:   | GA217553G1, Ccdc185 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217553G2, Ccdc185 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217553G3, Ccdc185 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | 4922505E12Rik                                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)