

Product datasheet for **GA217568**

Mouse Kprp activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kprp CRISPRa kit - CRISPR gene activation of mouse keratinocyte expressed, proline-rich
Symbol:	Kprp
Components:	GA217568G1, Kprp gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217568G2, Kprp gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217568G3, Kprp gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RUR4</u>
Synonyms:	1110001M24Rik; AA589586
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)