


Product datasheet for GA217575

Mouse Spag8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spag8 CRISPRa kit - CRISPR gene activation of mouse sperm associated antigen 8
Symbol:	Spag8
Components:	GA217575G1, Spag8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217575G2, Spag8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217575G3, Spag8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3V0Q6
Synonyms:	MH-SPAG8; SMP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)