

Product datasheet for **GA217684**

Mouse Tnni3k activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnni3k CRISPRa kit - CRISPR gene activation of mouse TNNI3 interacting kinase
Symbol:	Tnni3k
UniProt ID:	Q5GIG6
Synonyms:	Cark; D830019J24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)