

Product datasheet for **GA217715**

Mouse Xrra1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Xrra1 CRISPRa kit - CRISPR gene activation of mouse X-ray radiation resistance associated 1
Symbol:	Xrra1
Components:	GA217715G1, Xrra1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217715G2, Xrra1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217715G3, Xrra1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3U3V8</u>
Synonyms:	AI449753
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)