


## Product datasheet for **GA217763**

### Mouse Bdp1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Bdp1 CRISPRa kit - CRISPR gene activation of mouse B double prime 1, subunit of RNA polymerase III transcription initiation factor IIIB                                                                                                                                                                 |
| Symbol:       | Bdp1                                                                                                                                                                                                                                                                                                    |
| Components:   | GA217763G1, Bdp1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217763G2, Bdp1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217763G3, Bdp1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q571C7</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | AI662220; AW049390; AW456745; B130055N23Rik; G630013P12Rik; mKIAA1241; TAF3B1; TFC5                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)