

## Product datasheet for **GA217766**

### Mouse Cma2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Cma2 CRISPRa kit - CRISPR gene activation of mouse chymase 2, mast cell                                                                                                                                                                                                                                 |
| Symbol:       | Cma2                                                                                                                                                                                                                                                                                                    |
| Components:   | GA217766G1, Cma2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217766G2, Cma2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217766G3, Cma2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | Mcp10                                                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)