


## Product datasheet for **GA217772**

### Mouse Baiap3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Baiap3 CRISPRa kit - CRISPR gene activation of mouse BAI1-associated protein 3
Symbol:	Baiap3
Components:	GA217772G1, Baiap3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217772G2, Baiap3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217772G3, Baiap3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Bap3; Gm937
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)