

Product datasheet for **GA217800**

Mouse Shisal2a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Shisal2a CRISPRa kit - CRISPR gene activation of mouse shisa like 2A
Symbol:	Shisal2a
Components:	GA217800G1, Shisal2a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217800G2, Shisal2a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217800G3, Shisal2a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2A9G7
Synonyms:	OTTMUSG00000008243
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)