

Product datasheet for **GA217842**

Mouse Tmem235 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem235 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 235
Symbol:	Tmem235
UniProt ID:	B1AQL3
Synonyms:	Cldn27; Gm12581; OTTMUSG00000007589
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)