

## Product datasheet for **GA217857**

### Mouse H2alle activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | H2alle CRISPRa kit - CRISPR gene activation of mouse H2A histone family member LIE                                                                                                                                                                                                                            |
| <b>Symbol:</b>       | H2alle                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA217857G1, H2alle gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217857G2, H2alle gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217857G3, H2alle gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | Gm14484; OTTMUSG00000016790                                                                                                                                                                                                                                                                                   |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

