


Product datasheet for **GA217940**

Mouse Siglec15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Siglec15 CRISPRa kit - CRISPR gene activation of mouse sialic acid binding Ig-like lectin 15
Symbol:	Siglec15
Synonyms:	Cd33l3; EG620235; SIGLEC-I
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)