

Product datasheet for GA217969

Mouse Tmem170b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem170b CRISPRa kit – CRISPR gene activation of mouse transmembrane protein 170B
Symbol:	Tmem170b
Components:	GA217969G1, Tmem170b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217969G2, Tmem170b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217969G3, Tmem170b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P86050</u>
Synonyms:	EG621976
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

