

Product datasheet for **GA217983**

Mouse Kctd21 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kctd21 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 21
Symbol:	Kctd21
Components:	GA217983G1, Kctd21 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217983G2, Kctd21 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217983G3, Kctd21 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3URF8</u>
Synonyms:	EG622320
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)