


Product datasheet for **GA218005**

Mouse Dok6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dok6 CRISPRa kit - CRISPR gene activation of mouse docking protein 6
Symbol:	Dok6
Components:	GA218005G1, Dok6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218005G2, Dok6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218005G3, Dok6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q2MHE5
Synonyms:	Dok-6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)