

## Product datasheet for **GA218044**

### Mouse Gm6583 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Gm6583 CRISPRa kit - CRISPR gene activation of mouse predicted gene 6583                                                                                                                                                                                                                                      |
| <b>Symbol:</b>       | Gm6583                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA218044G1, Gm6583 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA218044G2, Gm6583 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA218044G3, Gm6583 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>1 CRISPRa-Enhancer vector, SKU GEI00056<br>1 CRISPRa scramble vector, SKU GEI00077 |
| <b>Synonyms:</b>     | EG625424                                                                                                                                                                                                                                                                                                      |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

