

Product datasheet for **GA218155**

Mouse Spint3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spint3 CRISPRa kit - CRISPR gene activation of mouse serine peptidase inhibitor, Kunitz type, 3
Symbol:	Spint3
Synonyms:	Gm14315; OTTMUSG00000016288
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)