

## Product datasheet for **GA218194**

### Mouse Zfp808 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Zfp808 CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 808                                                                                                                                                                                                                                  |
| <b>Symbol:</b>       | Zfp808                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA218194G1, Zfp808 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA218194G2, Zfp808 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA218194G3, Zfp808 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>1 CRISPRa-Enhancer vector, SKU GEI00056<br>1 CRISPRa scramble vector, SKU GEI00077 |
| <b>Synonyms:</b>     | EG630579; Gm7036                                                                                                                                                                                                                                                                                              |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

