

Product datasheet for GA218312

Mouse Tmem238 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem238 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 238
Symbol:	Tmem238
Components:	GA218312G1, Tmem238 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218312G2, Tmem238 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218312G3, Tmem238 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A9JSM3
Synonyms:	2210411K11Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

