


Product datasheet for **GA218340**

Mouse Zscan4f activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zscan4f CRISPRa kit - CRISPR gene activation of mouse zinc finger and SCAN domain containing 4F
Symbol:	Zscan4f
Components:	GA218340G1, Zscan4f gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218340G2, Zscan4f gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218340G3, Zscan4f gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3URS2</u>
Synonyms:	C330039H09Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)