

Product datasheet for GA218347

Mouse Cyp4a3l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp4a3l CRISPRa kit – CRISPR gene activation of mouse cytochrome P450, family 4, subfamily a, polypeptide 3l
Symbol:	Cyp4a3l
Components:	GA218347G1, Cyp4a3l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218347G2, Cyp4a3l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218347G3, Cyp4a3l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	RP23-118K16.4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

