

## Product datasheet for **GA218365**

### Mouse Ms4a4a activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Ms4a4a CRISPRa kit - CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 4A |
| Symbol:       | Ms4a4a   |
| Synonyms:     | EG666907   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                  |

#### Product images:



[View online >](#)