

## Product datasheet for **GA218379**

### Mouse *Ifit1b1* activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | <i>Ifit1b1</i> CRISPRa kit - CRISPR gene activation of mouse interferon induced protein with tetratricopeptide repeats 1B like 1 |
| Symbol:       | <i>Ifit1b1</i>   |
| Synonyms:     | Gm14446; <i>Ifit1c</i> ; OTTMUSG00000016644  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online >](#)