


Product datasheet for **GA218595**

Mouse Sap25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sap25 CRISPRa kit - CRISPR gene activation of mouse sin3 associated polypeptide
Symbol:	Sap25
Components:	GA218595G1, Sap25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218595G2, Sap25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218595G3, Sap25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q1EHW4</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)