

Product datasheet for **GA218710**

Mouse Hottip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hottip CRISPRa kit - CRISPR gene activation of mouse Hoxa distal transcript antisense RNA
Symbol:	Hottip
Components:	GA218710G1, Hottip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA218710G2, Hottip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA218710G3, Hottip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)