


Product datasheet for **GA219306**

Mouse Esp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Esp4 CRISPRa kit - CRISPR gene activation of mouse exocrine gland secreted peptide 4
Symbol:	Esp4
Components:	GA219306G1, Esp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA219306G2, Esp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA219306G3, Esp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm20580
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)