

Product datasheet for GA219782

Mouse B230354K17Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	B230354K17Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA B230354K17 gene
Symbol:	B230354K17Rik
Components:	GA219782G1, B230354K17Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA219782G2, B230354K17Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA219782G3, B230354K17Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)