

Product datasheet for **GA219788**

Mouse Lincmd1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lincmd1 CRISPRa kit - CRISPR gene activation of mouse long intergenic non-protein coding RNA of muscle differentiation 1
Symbol:	Lincmd1
Components:	GA219788G1, Lincmd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA219788G2, Lincmd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA219788G3, Lincmd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)