

Oligo annealing buffer (10X)

Catalog No. GE100007
Size 0.2 ml

Description
Oligo annealing buffer.

Appearance
Clear.

Form
Liquid, 10X Concentrate.

Shipping
Room temperature

Storage
Store at 4°C.

Stability
Stable for a minimum of one year from date of receipt when stored at 4°C.
1X working solution is stable for up to 4 weeks when stored at 4°C.

Related Products	Catalog No.
pCas-Guide cloning kit	GE100001
pCas-Guide plasmid	GE100002
pCas-Scramble	GE100003
CF3, the forward sequencing primer	GE100008

Complete genome editing tools OriGene offered, www.origene.com/cas9

For Research Use Only. Not for use in diagnostic procedures.