

CF3, gRNA sequencing primer

Catalog No.
GE100008

Size
100 pmoles

Description

CF3, the forward sequencing primer for CRISPR vectors to sequence target sequence cloned; it can be used for all gRNA CRISPR vectors except pT7-Guide-IVT, including all-on-one and gRNA only cloning.

Primer sequence: 5'-ACGATACAAGGCTGTTAGAGAG-3'

Appearance

Lyophilized. Please resuspend in 10 µL H₂O to make a 10 µM solution.

Shipping

Room temperature

Storage

Store at 4°C.

Stability

Stable for a minimum of one year from date of receipt when stored at 4°C.
Reconstituted solution is stable for up to 4 weeks when stored at 4°C.

Related Products

pCas-Guide cloning kit

pCas-Guide plasmid

pCas-Scramble

Oligo annealing buffer

Complete genome editing tools OriGene offered, www.origene.com/cas9

Catalog No.

GE100001

GE100002

GE100003

GE100007

For Research Use Only. Not for use in diagnostic procedures.