

## Product datasheet for **GE100010**

### pLenti-Cas-Guide

#### Product data:


Product Type: CRISPR Vectors

Function: All-in-One CRISPR/Cas9 vectors

Features:

- Lentiviral vector, ideal for hard-to-transfect and animal models
- All-in-one vector
- 1. Expresses human codon-optimized Cas9, C-terminal Myc-DDK tagged
- 2. BamH I and BsmB I can be used for genomic target sequence cloning
- 3. gRNA expression driven by U6 promoter
- 4. Lenti particles can be produced using the high efficient [lenti packaging kit](#)

#### Product images:



[View online »](#)