

Product datasheet for **GE100014**


pT7-Cas9

Product data:

- Product Type:** CRISPR Vectors
- Function:** CRISPR/Cas9 T7 Vectors
- Features:**
- T7 driven Cas9 expression vector
 - 1. Cas9 mRNA can be produced *via in vitro* transcription.
 - 2. For mRNA mediated genome editing via microinjection or mRNA transfection, both Cas9 mRNA and gRNA transcript are needed. T7 driven gRNA can be produced from a separate vector, SKU GE100025.

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)