

Product datasheet for **GE100020**

pT7-Cas9-Nickase (D10A)

Product data:

Product Type: CRISPR Vectors


Function: Cas9 Nickase

Features: T7 driven Cas9-Nickase expression vector

- Serves as a template for T7 *in vitro* transcription of Cas9 D10A mRNA
- Cas9 D10A expressed contains a C-terminal Myc-DDK tag
- T7 driven gRNA can be produced from a separate vector, GE100025

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)