

Product datasheet for **GE100022**

pCas-Guide-EF1a-CD4

Product data:

Product Type: CRISPR Vectors


Function: All-in-One CRISPR/Cas9 vectors

Features: All-in-one vector

- Expresses human codon-optimized Cas9, C-terminal Myc-DDK tagged
- BamH I and BsmB I sites are used for genomic target sequence cloning
- Expresses gRNA containing the inserted target sequence and the gRNA scaffold
- EF1a promoter driven CD4 to enrich transfected cells
- CF3 primer (GE100008) (GE100008) can be used to sequence the target sequence after cloning

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)