

Product datasheet for GE100029

pLenti-Cas9-IRES-Puro

Product data:

Product Type: CRISPR Vectors

Function: Cas9 only vectors

Features: Cas9 only vector with puro selection in lenti backbone

• 3rd generation Lenti vector, ideal for hard-to-transfect cells

• Cas9 under CMV promoter for mammalian expression after transfection or transduction

• IRES-Puro is within the Lentiviral packaging cassette and would allow Puromycin selection in


both the transfected and the transduced cells.

• Lenti particles can be produced using high efficiencylenti packaging kits

Disclaimer:

These products are manufactured and supplied by OriGene under license from ERS.

Product images:





OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com