

Product datasheet for **GE100031**


pLenti-Cas9-P2A-tGFP

Product data:

- Product Type:** CRISPR Vectors
- Function:** Cas9 only vectors
- Features:** Cas9 only vector with GFP expression in lenti backbone
- 3rd generation Lenti vector, ideal for hard-to-transfect cells
 - Cas9 under CMV promoter for mammalian expression after transfection or transduction
 - P2A-tGFP is within the Lentiviral packaging cassette and would allow GFP expression in both the transfected and the transduced cells.
 - Lenti particles can be produced using high efficiency [lenti packaging kits](#)

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)