

# Product datasheet for GE100032

### pLenti-Guide-Puro

#### **Product data:**

Product Type:

Function:

Features:

CRISPR Vectors gRNA only vectors

gRNA only vector in lenti backbone

- gRNA expression under U6 promoter
- 3rd generation Lenti vector, ideal for hard-to-transfect cells
- BamH1 and BsmB I sites are used for target sequence cloning
- Puro resistant gene allows puro selection in both transfected and transduced cells
- Lenti particles can be produced using high efficiency<u>lenti packaging kits</u>
- CF3 primer (GE100008) can be used to sequence the target sequence after cloning
- E. Coli selection: chloramphenicol (34 ug/mL)3rd generation Lenti vector, ideal for hard-totransfect cells

## **Product images:**





This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn