

Product datasheet for **GE100038**


Cas9 insertion kit at AAVS1 locus (Puro selection)

Product data:

Product Type:	CRISPR Vectors
Function:	CRISPR/AAVS1 kit
Features:	<p>Complete vector kit for targeted Cas9 insertion into human AAVS1 site via CRISPR</p> <p>Components:</p> <ul style="list-style-type: none"> • pCas-Guide-AAVS1 (GE100023) (SKU=GE100023) • pAAVS1-Cas9-Puro-DNR (GE100037) (SKU GE100037) • pCas-Guide-scramble (GE100003) (SKU=GE100003) <p>Cas9 driven by CMV promoter</p> <p>PGK-Puro allows mammalian puromycin stable selection</p>

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)