

Product datasheet for **GE100039**


pAAVS1-Cas9-BSD-DNR

Product data:

- Product Type:** CRISPR Vectors
- Function:** Cas9 only vectors
- Features:** Cas9 in AAVS1-BSD donor vector
- CMV-Cas9 and PGK-BSD in AAVS1 donor vector.
 - Contains AAVS1 left and right homologous sequences (550 bp each side)
 - CMV-Cas9 and PGK-BSD cassette can be integrated into human AAVS1 site when cotransfected with pCas-Guide-AAVS1 (GE100023)
 - When transfected into cells alone, can be used as regular Cas9 mammalian expression vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)