

## Product datasheet for **GE100040**


### Cas9 insertion kit at AAVS1 locus (BSD selection)

#### Product data:

Product Type:	CRISPR Vectors
Function:	CRISPR/AAVS1 kit
Features:	<p>Complete vector kit for targeted Cas9 insertion into human AAVS1 site via CRISPR</p> <p>Components:</p> <ul style="list-style-type: none"> <li>• pCas-Guide-AAVS1 (GE100023) (SKU GE100023)</li> <li>• pAAVS1-Cas9-BSD-DNR (GE100039) (SKU GE100039)</li> <li>• pCas-Guide-Scramble (GE100003) (SKU GE100003)</li> </ul> <p>Cas9 driven by CMV promoter</p> <p>PGK-BSD allows mammalian blasticidin stable selection</p>

**Disclaimer:** These products are manufactured and supplied by OriGene under license from ERS.

#### Product images:



[View online »](#)