

## Product datasheet for **GE100042**

### pGuide

#### Product data:

Product Type: CRISPR Vectors


Function: gRNA only vectors

Features: gRNA only vector (Cas9 not present)

- gRNA expression under U6 promoter
- BamH I and BsmB I sites are used for target sequence cloning
- Expresses gRNA containing the inserted target sequence and the gRNA scaffold
- CF3 primer (GE100008) (GE100008) can be used to sequence the target sequence after cloning

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

#### Product images:



[View online »](#)