

Product datasheet for **GE100044**


pGuide-EF1a-GFP

Product data:

- Product Type:** CRISPR Vectors
- Function:** gRNA only vectors
- Features:**
- gRNA only vector (Cas9 not present)
 - gRNA expression under U6 promoter
 - BamH I and BsmB I sites are used for target sequence cloning
 - Expresses gRNA containing the inserted target sequence and the gRNA scaffold
 - GFP as a reporter for tracking or sorting transfected cells
 - CF3 primer (GE100008) (GE100008) can be used to sequence the target sequence after cloning

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)