

## Product datasheet for **GE100045**


### pLenti-EF1a-Cas-Guide

#### Product data:

- Product Type:** CRISPR Vectors
- Function:** All-in-One CRISPR/Cas9 vectors
- Features:**
- Lentiviral vector, ideal for hard-to-transfect and animal models
  - All-in-one vector
1. Expresses human codon-optimized Cas9, C-terminal Myc-DDK tagged
  2. BamH I and BsmB I can be used for genomic target sequence cloning
  3. gRNA expression driven by U6 promoter
  4. Lenti particles can be produced using the high efficient [lenti packaging kit](#)

**Disclaimer:** These products are manufactured and supplied by OriGene under license from ERS.

#### Product images:



[View online »](#)