

Product datasheet for **GE100046**

pAAVS1-EF1a-puro-DNR CRISPR Vector


Product data:

Product Type: CRISPR Vectors

Function: AAVS1 donor vector

- Features:
- Donor vector for targeted transgene insertion into human AAVS1 site
 - Contains left and right AAVS1 homologous sequences (550 bp each side)
 - Multiple cloning sites for transgene cloning, compatible for shuttling from the TrueORF clones
 - The expression of the target gene will be under EF1a promoter, a constitutive promoter as an alternative of CMV promoter. the latter could be silenced in certain cellular environment
 - Puromycin mammalian selection under PGK promoter
 - Targeted insertion when cotransfected with pCas-Guide-AAVS1 (SKU GE100023)

Product images:



[View online »](#)

