

Product datasheet for **GE100049**

AAVS1 Transgene Knockin kit (EF1a, BSD)

Product data:

- Product Type:** CRISPR Vectors
- Function:** CRISPR/AAVS1 kit
- Features:** Complete vector kit for targeted transgene insertion into human AAVS1 site via CRISPR
Components:
- pCas-Guide-AAVS1 (SKU GE100023)
 - pAAVS1-EF1a-BSD-DNR (SKU GE100048)
 - pCas-Guide-scramble (SKU GE100003)
- Transgene needs to be cloned into pAAVS1-EF1a-BSD-DNR

Product images:

