


Product datasheet for **GE100051**

pROSA26-Puro-DNR

Product data:

- Product Type:** CRISPR Vectors
- Function:** ROSA26 donor vector
- Features:**
- Donor vector for targeted transgene insertion into Mouse ROSA26 locus
 - Contains left and right ROSA26 homologous sequences (550 bp each side)
 - Multiple cloning sites for transgene cloning, compatible for shuttling from the TrueORF clones
 - Puromycin mammalian selection under PGK promoter
 - Targeted insertion when cotransfect with pCas-Guide-ROSA26 (GE100050) (SKU GE100050)

Product images:



[View online »](#)

EcoRI
BamHI KpnI
RBS
Kozak
Consensus
SgfI

CTATAGGGCGGCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

AscI
MluI
NotI

GGCGCGCCAGATCTCAAGCTTAACTAGTTAGCGGACCGACGCGTACGCGG CCG CTC GAG CAG
 T R P L E Q

Myc. Tag
DDK. Tag

AAA CTC ATC TCA GAA GAG GAT CTG GCA GCA AAT GAT ATC CTG GAT TAC AAG GAT GAC GAC
 K L I S E E D L A A N D I L D Y K D D D

PmeI

GAT AAG GTT TAA AC
 D K V Stop