

Cas9 In Vitro Digestion Buffer (10X)

Catalog No.

GE100053

Size

0.2 ml

Description

Cas9 in vitro digestion buffer, 10x

Appearance

Clear

Form

Liquid, 10X Concentrate

Shipping

Room temperature

Storage

Store at -20°C

Stability

Stable for a minimum of one year from date of receipt when stored at -20°C.
1X working solution is stable for up to 4 weeks when stored at 4°C.

Related Products

Cas9 purified protein

Catalog No.

TP790148

For Research Use Only. Not for use in diagnostic procedures.