


Product datasheet for **GE100054**

Cas9 mRNA

Product data:

Product Type:	CRISPR Accessories
Description:	Purified Cas9 mRNA (<i>Streptococcus Pyogenes</i>), capped and polyA-tailed, 10 ug
Recommended Application:	Transfection or microinjection with gRNA for genome editing
Species:	<i>Streptococcus Pyogenes</i>
Features:	Capped and polyA-tailed

Product images:



Cas9 mRNA with or without sgRNA were transfected into HEK 293T cells using Viromer mRNA. After 48 hrs, cells were harvested and genome editing was measured using T7E1 assay.



[View online >](#)